



NANO DSC NANO DIFFERENTIAL SCANNING CALORIMETER



Getting Started Guide Models 602000 and 602001

Notice

The material contained in this manual, and in the online help for the software used to support this instrument, is believed adequate for the intended use of the instrument. If the instrument or procedures are used for purposes other than those specified herein, confirmation of their suitability must be obtained from TA Instruments. Otherwise, TA Instruments does not guarantee any results and assumes no obligation or liability. TA Instruments also reserves the right to revise this document and to make changes without notice.

TA Instruments may have patents, patent applications, trademarks, copyrights, or other intellectual property covering subject matter in this document. Except as expressly provided in written license agreement from TA Instruments, the furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property.

TA Instruments Operating Software, as well as Module, Data Analysis, and Utility Software and their associated manuals and online help, are proprietary and copyrighted by TA Instruments. Purchasers are granted a license to use these software programs on the module and controller with which they were purchased. These programs may not be duplicated by the purchaser without the prior written consent of TA Instruments. Each licensed program shall remain the exclusive property of TA Instruments, and no rights or licenses are granted to the purchaser other than as specified above.

TA Instruments can accept no liability for loss or damage, however caused, arising from the faulty or incorrect use of its products. TA Instruments shall not be liable for any damages caused by interactions between exogenous materials (e.g. chemicals) and parts of the instrument. This includes interactions of gaseous, liquid or solid materials with, for instance, ampoule surfaces and/or parts of the instrument. It also includes gases or vapors leaking from ampoules (e.g. originating from chemical reactions producing gaseous substances), with subsequent cause of damage to the calorimeter.

©2013-2017 by
TA Instruments — Waters LLC
159 Lukens Drive
New Castle, DE 19720

Introduction

Important: TA Instruments Manual Supplement

Please click the [TA Manual Supplement](#) link to access the following important information supplemental to this Getting Started Guide:

- TA Instruments Trademarks
- TA Instruments Patents
- Other Trademarks
- TA Instruments End-User License Agreement
- TA Instruments Offices

Notes, Cautions, and Warnings

This manual uses NOTES, CAUTIONS, and WARNINGS to emphasize important and critical instructions. In the body of the manual these may be found in the shaded box on the outside of the page.

NOTE: A NOTE highlights important information about equipment or procedures.

CAUTION: A CAUTION emphasizes a procedure that may damage equipment or cause loss of data if not followed correctly.

MISE EN GARDE: UNE MISE EN GARDE met l'accent sur une procédure susceptible d'endommager l'équipement ou de causer la perte des données si elle n'est pas correctement suivie.

A WARNING indicates a procedure that may be hazardous to the operator or to the environment if not followed correctly.

Un AVERTISSEMENT indique une procédure qui peut être dangereuse pour l'opérateur ou l'environnement si elle n'est pas correctement suivie.

Regulatory Compliance

Safety Standards

EMC Directive

This instrument has been tested to meet the European Electromagnetic Compatibility Directive (EMC Directive, 2014/30/EU). The Declaration of Conformity for your instrument lists the specific standards to which the unit was tested.

The instrument was designed specifically as a test and measuring device. Compliance to the EMC directive is through IEC 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements (1998).

As noted in the IEC 61326-1, the instrument can have varying configurations. Emissions may, in non-typical applications, exceed the levels required by the standard. It is not practical to test all configurations, as the manufacturer has no control over the user application of the instrument.

Immunity Testing

The instrument was tested to the requirements for laboratory locations.

Emission Testing

The instrument fulfills the limit requirements for Class A equipment but does not fulfill the limit requirements for Class B equipment. The instrument was not designated to be used in domestic establishments.

Low Voltage Directive (Safety)

In order to comply with the European Low Voltage Directive (2014/35/EU), this equipment has been designed to meet IEC 1010-1 (EN 61010-1) standards. To comply with requirements in the USA, this instrument has been tested to the requirements of UL61010a-1.

Safety

Electrical Safety

There are no user-serviceable parts inside the instrument; voltages as high as 125/250 VAC are present in this system.

Lifting the Instrument

The Nano DSC is not a portable instrument. In order to avoid injury, particularly to the back, please follow this advice:

WARNING: Use appropriate care when unpacking or moving the instrument. It may be too heavy for some individuals working alone to handle safely.

AVERTISSEMENT: Soyez prudent lors du dépotage ou du déplacement de l'instrument. Il peut être trop lourd à manipuler en toute sécurité pour des personnes travaillant seules.

DANGER: High voltages are present in this instrument. Maintenance and repair of internal parts must be performed only by TA Instruments or other qualified service personnel.

DANGER: Présence de tensions élevées dans cet instrument. La maintenance et la réparation des pièces internes doivent être effectuées uniquement par TA Instruments ou tout autre personnel d'entretien qualifié.

Table of Contents

Introduction	3
Important: TA Instruments Manual Supplement	3
Notes, Cautions, and Warnings	4
Regulatory Compliance	4
Safety Standards	4
Safety	5
Electrical Safety	5
Lifting the Instrument	5
Chapter 1: Introducing the Nano DSC	8
Overview	8
Instrument Models Covered in this Guide	8
Applications	9
System Components	10
Options and Accessories	10
Instrument Specifications	11
Chapter 2: Installing the Nano DSC	12
Unpacking/Repacking the Nano DSC	12
Installing the Instrument	12
Inspecting the System	12
Choosing a Location	13
In	13
On	13
Near	13
Away from	13
Power Requirements	13
Gas Requirements	14
Electrical Connections	14
Starting the Nano DSC	15
Shutting Down the Nano DSC	16
Chapter 3: Use, Maintenance, & Diagnostics	17
Overview	17
Preparing the Buffer Solutions	17
Preparing Dilute Biopolymer Solutions	18
Degassing the Buffer and Sample Solutions	18
Running Baseline and Sample Solution Scans	19
Selecting Solutions	19
Filling Cylindrical Cells	19

Filling the Capillary Cells	20
Cell Pressurization and Checking for Bubbles	23
Cleaning the Nano DSC	25
Cylindrical Cell Design Flushing Equipment	26
Capillary Cell Design Flushing Equipment	27
Software	28
Maintaining the Nano DSC	29
Balance and Residual	29
Purging the Nano DSC	31
Cleaning O-Rings	31
Cleaning the Instrument	31
Regular Cleaning	31
Thorough Cleaning	31
Troubleshooting the Instrument	33
Appendix A: Pressure Perturbation Experiments	35
Introduction	35
PPC with DSCRun	35
Experimental Set-up	37
PPC Fitting with NanoAnalyze	37
References	40
Appendix B: Guardian Option	41
Guardian Option	41

Chapter: 1

Introducing the Nano DSC

Overview

The TA Instruments Nano Differential Scanning Calorimeter (Nano DSC) is designed specifically for the measurement of absolute heat capacities of biopolymers in dilute solution.




Figure 1 The Nano DSC instrument.

The Nano DSC offers improvement in sensitivity and baseline repeatability. The instrument is a power compensation design using a completely solid-state thermostat and is equipped with fixed continuous platinum capillary or 24K gold cylindrical cells.

With an operating temperature range of -10 to 130°C (or 160°C for high-temperature units), scan rates up to $2^{\circ}\text{C}/\text{min}$ in both the heating and cooling directions, and baseline repeatability of $\pm 0.4 \mu\text{cal}/\text{C}$ @ $1^{\circ}\text{C}/\text{min}$. The Nano DSC can be used to study the thermal denaturation of many proteins with as little as $52 \mu\text{g}$ or less of sample.

Instrument Models Covered in this Guide

Models 602000 and 602001

Applications

The Nano DSC allows you to study almost any thermal transition occurring in dilute biopolymer solution in the temperature range of -10 to 130°C . It also offers a new dimension in sensitivity, baseline noise, and baseline repeatability in both heating and cooling modes.

The measurement of absolute heat capacities using the instrument enables the examination of biopolymer solution conformation and solvation, as well as subtle changes in folding.

Studies of biopolymer thermal stability and the effects of biopolymer or ligand interactions on macromolecular structure can be accomplished with the Nano DSC. The baseline repeatability ensures that reliable DCp data can be obtained from a single scan.

Typical uses for the Nano DSC include:

- Biopolymer solution conformation and solvation (absolute heat capacities)
- Biopolymer stability (protein denaturation)
- Biopolymer structure (domain organization)
- Bioengineering (mutant proteins)
- Ligand interactions (drug binding to proteins or nucleic acids)
- Membrane structure (lipid bilayers, membrane proteins)
- Polynucleotides (helix to coil transitions)

System Components

The following items make up a Nano DSC system:

- Nano Differential Scanning Calorimeter
- Power cord
- Computer system (optionally available from TA Instruments) and display monitor with accessories, and all required software with online Help and manual. The computer system is current technology at the time of purchase, and uses the Microsoft Windows operating system.
- Pressure handle and 1 replacement o-ring.
- 1-mL filling syringe and filling needles (cylindrical cell instrument only)
- Cell cleaning adapter with silicone tubing (cylindrical cell instrument only)
- 2.5-mL micropipette with tips and silicone tubing (capillary cell instrument only)
- Silicone tubing for cleaning (capillary cell instrument only)
- Windows-based Experiment Setup, Data Collection and Analysis Software (This software was previously installed at the factory and the supplied disk is typically used for software reinstallation.)

Options and Accessories

The following items are available to complement your Nano DSC instrument:

- Vacuum Degassing System
- Personal computer with Microsoft Windows Operating System
- Additional filling needles
- Replacement 1-mL syringe
- Micropipetter and micropipette tips
- Replacement Nano DSC cleaning adapter
- Replacement silicone tubing for cleaning adapter
- Replacement pressure handle
- Replacement pressure o-ring kit

Instrument Specifications

The table found below contains the technical specifications for the Nano DSC instrument.

Table 1: Nano DSC Technical Specifications

Item/Area	Specifications
Dimensions	Depth 53 cm (21 in.) Width 35 cm (14 in.) Height 28 cm (11 in.)
Weight	17 kg (37 lbs)
Power	100–240 VAC, 3 amps. 50 or 60 Hz
Electrical power cord	The plug of the cord must be rated to carry at least 125% of the product current rating. The cord length must be less than 4.5 meters and must be UL or CSA approved.
Operating environmental conditions	Temperature: 15°C to 30°C Relative humidity: 5 to 80% (non-condensing) Installation category II Pollution degree 2 Maximum altitude: 2000 m (6560 ft)
Emissions class	Class A
Temperature range	–10 to 130°C (Standard Temperature) or –10 to 160°C (High Temperature)
Temperature scan rate	0°C to 2°C/min
Response time	Halftime = 5 s
Baseline repeatability	0.028 μ W
Baseline noise	Standard deviation 0.015 μ W
Cell construction	Capillary in platinum or cylindrical in 24K gold
Cell volume	0.30 mL nominal capillary or 0.33 mL cylindrical
Cell pressurization range	0 to 6 atm

Chapter: 2

Installing the Nano DSC

Unpacking/Repacking the Nano DSC

WARNING: Have an assistant help you unpack this unit. Do not attempt to do this alone.

AVERTISSEMENT: Faites-vous aider par une personne pour dépoter cet appareil. N'essayez pas de le faire tout seul.

CAUTION: To avoid mistakes, read this entire chapter before you begin installation.

MISE EN GARDE: Pour éviter de commettre des erreurs, lisez tout le chapitre avant de commencer l'installation.

The instructions needed to unpack and repack the instrument are found as separate unpacking instructions in the shipping box. You may wish to retain all of the shipping hardware and boxes from the instrument in the event you wish to repack and ship your instrument.

Installing the Instrument

Before shipment, the instrument is inspected both electrically and mechanically so that it is ready for operation upon proper installation. Only limited instructions are given in this manual; consult the online documentation for additional information. Installation involves the following procedures:

- Inspecting the system for shipping damage and missing parts
- Connecting the Nano DSC to the computer
- Connecting USB cables

It is recommended that you have your Nano DSC installed by a TA Instruments Service Representative; call for an installation appointment when you receive your instrument.

Inspecting the System

When you receive your instrument, look over the instrument and shipping container carefully for signs of shipping damage, and check the parts received against the enclosed shipping list.

- If the instrument is damaged, notify the carrier and TA Instruments immediately.
- If the instrument is intact but parts are missing, contact TA Instruments.

Choosing a Location

It is important to choose a location for the instrument using the following guidelines. The Nano DSC should be:

In

- a temperature- and humidity-controlled area. Temperatures should be in range 15 to 30°C.
- a clean, vibration-free environment, preferably on the ground floor in the building. It should be located away from pumps, motors, or other devices which produce vibrations.
- an area with ample working and ventilation space. At least 18 by 18 inches is needed for the instrument. Additional space is needed for the computer and (if present) printer.

On

- a stable work surface.

Near

- a power outlet. See the [“Power Requirements”](#) section below.
- your computer.

Away from

- dusty environments.
- exposure to direct sunlight.
- direct air drafts (fans, room air ducts).
- poorly ventilated areas.
- noisy or mechanical vibrations.
- high traffic areas, where constant movements from passing personnel could create air currents or mechanical disturbances.

Power Requirements

The Nano DSC requires a grounded, single-phase power source. A three-conductor line cord ensures a safety ground. The operating voltage and line frequency were preset at the factory for 100–240 VAC, 50 or 60 Hz operation.

It is highly recommended that the instrument and computer system should be operated from a single power source. Use a power distribution strip with a surge suppressor function. If the local power mains experience fluctuating voltage or surge conditions, it is also advisable to use a power conditioning uninterruptible power supply (UPS).

The instrument requires one (1) 15-A line for the instrument and computer system. Attach all items to a single surge suppressor plug strip. Plug the surge suppressor strip into an isolated power line. An isolated power line is one that is used only for electronic instruments (i.e., no inductive motors, compressors, or heaters) that do not create power surges upon start-up.

Attach the instrument and all computer accessories to the surge suppressor.

Gas Requirements

Because the Nano DSC is capable of running below room temperature, the calorimeter's thermal shield features a sealed environment that has been filled with dry nitrogen gas.

CAUTION: The purge port valve on the back of the Nano DSC should remain in the closed position at all times to maintain the integrity of the nitrogen purge.

MISE EN GARDE: La vanne de l'orifice de drainage située à l'arrière du Nano DSC [ITC] doit toujours rester en position fermée pour maintenir l'intégrité de l'azote drainé.

NOTE: Purging of the canister is not a routine maintenance operation; contact TA Instruments for service if the Nano DSC exhibits unstable baselines especially above 50 or 60°C.

Electrical Connections

- 1 Make sure all of the equipment is turned off.
- 2 Make sure the power switch to the surge suppressor is in the OFF position.
- 3 Plug the power cord into the Nano DSC and the surge suppressor strip.
- 4 Assemble the computer system per the instructions included in its shipment carton.
- 5 Install DSCRun software, following the instructions in the “Installing DSC Run” document found on the TA Instruments software CD.
- 6 Connect the calorimeter to the computer system with the supplied USB cable.
- 7 Ensure that airflow to the cooling fan on the rear panel of the calorimeter is not blocked by bench back-splashes, walls, or other equipment.

Starting the Nano DSC

Once you have completely set up the calorimeter and computer system, start the instrument as follows:

- 1 Turn on the surge suppressor power switch and the computer system and monitor. Allow the computer to boot up.
- 2 Turn on the power switch to the calorimeter, which is located on the back panel. The front LED will light up when in the “on” position.



Figure 2 Power LED on Nano DSC instrument.

- 3 Start DSCRun™ software on the computer (see the figure below). Instrument control and regulation is supplied by the software. Initial thermal stabilization may take an hour or longer.




Figure 3 DSCRun software screen.

Shutting Down the Nano DSC

You can leave the instrument and its associated components on when the Nano DSC will be inactive for two to five days.

If the Nano DSC will be inactive for more than five days, TA Instruments recommends that you turn all equipment off.

Chapter: 3

Use, Maintenance, & Diagnostics

Overview

A typical Nano DSC experiment involves the following:

- Preparing the buffer/solvent and biopolymer solutions
- Running the baseline and the sample solution scans
- Cleaning the calorimeter
- Performing a thermodynamic analysis

Each step is briefly described below. Additional information is provided in the online Help supplied with the software program. It is assumed that you are familiar with standard laboratory procedures and techniques. It is critical that the Nano DSC cells be cleaned immediately at the end of each experiment. Although the calorimeter can be left idle with water-filled cells at 25°C when not performing experiments, it will generally perform better when it is left continuously scanning. This is done by filling the cells with deionized degassed water, pressurizing the cell to 45 psi (3 atm), and scanning up and down continuously between the upper and lower temperature limits at 1°C/min.

Preparing the Buffer Solutions

The solvent systems used for preparation of dilute biopolymer solutions are typically buffered aqueous solutions, which often contain supporting strong electrolytes (for example, NaCl or KCl) to adjust the ionic strength. To prepare a buffer, follow these steps:

- 1 Weigh an appropriate amount of the acidic or basic form of a weak acid (or weak base) into a measured volume of distilled or deionized water.
- 2 Adjust the pH of the buffer to the desired value by adding either strong acid (for example, HCl) or strong base (for example, NaOH), while monitoring the buffer pH using a pH electrode and meter. The buffer will provide the calorimeter baseline. It will also be used to prepare and dialyze the protein (or nucleic acid) solution.
- 3 Prepare a rather large volume of the buffer (for example, 1 L). It should be noted that buffers used in DSC experiments should be carefully chosen to meet the following criteria:
 - The pKa should be as independent as possible of the temperature.
 - The ΔH for proton ionization from the buffer acid should be small.
 - All components of the buffer/solvent solution should be thermally stable. (The buffer should not precipitate or change color on boiling.)

Preparing Dilute Biopolymer Solutions

Choose the concentration of the biopolymer to allow optimization of the amount of thermodynamic information. In the case of a typical protein, good heat capacity data can usually be measured for solutions with protein concentrations in the 0.2 to 5 mg/mL range. When possible, obtain the information regarding the optimum concentration and solution conditions for DSC studies of specific biopolymers from the literature.

The minimum of volume of dilute protein solution required to fill the Nano DSC cylindrical cell is 0.5 mL. About 0.65 mL is required to fill a capillary cell. A protein solution is normally prepared by weight or by dilution of a concentrated stock solution. For the best results, dialyze the protein (or other biopolymer sample) before your DSC experiment. Then use the dialysate for the reference solution.

It must be noted that high sensitivity DSC studies should be done with highly purified biopolymers. Accepted techniques should be used for sample purification prior to solution preparation and dialysis. It is critical to know the solute (biopolymer) molecular weight, concentration and state of oligomerization prior to attempting a thermodynamic analysis of the heat capacity data.

Degassing the Buffer and Sample Solutions

Typically, if a solution is heated, gas bubbles will form as the solubility of dissolved gases (such as O₂ and N₂) is decreased with increasing temperature. If gas bubble formation occurs in the DSC cells during the run, the resulting heat capacity data will be rather noisy since abrupt changes in the apparent heat capacity will result from the bubble driven liquid displacement effects.

All reference and sample solutions must be degassed prior to being placed in the DSC to minimize the possibility of gas bubble formation during the run. Pull a vacuum of 15 to 25 inches Hg on the solutions for a period of 10 to 15 min to degas a sample. An accessory degassing station is available from TA Instruments.

Gas bubble formation during the run is retarded by the cell pressure (up to 3 to 6 atm) provided by the DSC manostat. However, it is still necessary to remove the atmospheric gasses from the solutions.

Running Baseline and Sample Solution Scans

The second major step in an experiment includes running the baseline and sample scans. The first step in this procedure is to fill the cells. There are two protocols for filling the cells, depending on the type of Nano DSC cell design. Cylindrical cells utilize a syringe while capillary cells use a pipette. Follow the instructions that relate to your particular type of cell.

Selecting Solutions

Follow these guidelines when choosing the solution to use in your experiments:

- Baseline Run: Fill both the sample cell and reference cells with the same solution (buffer dialysate).
- Sample Run: Fill the sample cell with the dialyzed protein (or other biopolymer) solution and the reference cell with the dialysate.

Filling Cylindrical Cells

Filling the cell is one of the most important parts of conducting a DSC experiment and should be done with extreme care. The objective is to fill the cell with the appropriate solution without introducing any air bubbles. Even the smallest bubble can result in heat capacity data that is erratic and offset. You will need a 1-mL Hamilton syringe equipped with a 20-gauge needle (this needle and syringe are supplied with your instrument).

NOTE: Proteins in solution have a tendency to foam when manipulated with a pipette.

- 1 Make sure that the cells are empty. Starting with the reference cell (left-most access tube as viewed from the front of the DSC), begin by filling the syringe with approximately 0.5 mL of solution. You may find that larger volumes make the procedure easier, but it is possible to load the cell correctly with the lesser amount. Be careful not to trap any air bubbles.
- 2 Insert the syringe into the reference cell and load the calorimeter in the following manner ():
 - a Position the needle so that it **lightly** touches the bottom of the reference cell (left access tube) initially. See the figure below.



Figure 4 Inserting the syringe.

- b** Raise the syringe off the bottom approximately 1/16 in.
 - c** Gently inject the solution into the cell until you see it coming out of the access tube.
- 3** Slowly pump about 0.1 mL of solution in and out of the access tube and watch to see if bubbles are coming out of the access tube. Perform this operation a number of times until it is apparent that there are **no** bubbles in the cell.
- 4** Slowly withdraw the needle from the access tube while simultaneously injecting solution to clear away any bubbles trapped by the needle.
- 5** When you are finished, the reservoir located at the top of the access tube should be filled up to halfway with solution.
- 6** Repeat the entire process to fill the sample cell (right access tube). In the case of a baseline run, the sample cell will be filled with the same solution (buffer dialysate) that was placed in the reference cell. In the case of the sample run, the sample cell will be filled with the dialyzed protein (or other biopolymer) solution.

Filling the Capillary Cells

Equipment Required:

- Adjustable 500 to 2500 μL pipette (Eppendorf Series 2000 supplied)
- Tubing (1/16-in. ID Manosil silicone rubber)

NOTE: Additional care must be taken with proteins and other biopolymers since they have a tendency to foam when manipulated with a syringe.

As with a cylindrical cell, filling the capillary cell is one of the most important parts of conducting a DSC experiment and should be done with extreme care. The objective is to fill the cell with the appropriate solution without introducing any air bubbles. Even the smallest bubble can result in heat capacity data that is erratic and offset. You will need an adjustable volume pipette with the volume set between 800 and 1000 μL . (With experience, this volume can be reduced to 0.65 mL).

You will also need four clean pipette tips and four short lengths (approximately 0.5 in) of 1/16 in. I.D. Manosil silicone rubber tubing.

- 1** Be sure that the samples have been properly degassed. See [“Degassing the Buffer and Sample Solutions” on page 18.](#)
- 2** Displace any residual liquid from the cells by blowing them out with an inert gas (N_2) at low pressure.

- 3 Prepare the pipette tips, used with the micropipette for loading and cleaning the DSC cells. Fit the disposable plastic tips with a short (approximate 0.5-in. long) length of Manosil silicone tubing to interface with the cell access tubes. The figures below illustrate how to trim and prepare the pipette tips. A length of the appropriate silicone tubing is provided with each calorimeter. Additional tubing can be ordered through your laboratory supplier or TA Instruments.




Figure 5 Cutting the silicon tubing on the pipette tip (top image); completed pipette tip (bottom image).

- 4 For best bubble-free filling, leave the cells full of water while the instrument is idle. Empty the cells immediately before refilling. Aspirate by hand with a syringe, or with a vacuum pump just long enough to remove the liquid (about 2 seconds). This promotes wetting by preventing the cell surface from becoming completely dry.

Aspirate by hand with a syringe or with a vacuum pump just long enough to remove the liquid (about 2 seconds). This promotes wetting by preventing the cell surface from becoming completely dry.
- 5 Start with the reference cell (left-most pair of access tubes as viewed from the front of the DSC). Attach one pipette tip to either the front or rear access tube.
- 6 Attach the second pipette tip to the pipette and fill it with at least 650 μL of the reference solution (buffer dialysate).

- 7 Attach the silicone extension of the pipette tip to the reference cell's inlet access tube. See the below. Fill the cell slowly (over 5 to 10 seconds) by gently depressing the pipette plunger.



Figure 6 Pipette and silicone tubing attachment to reference cell access tubes.

- 8 Pipette and silicone tubing attachment to reference cell access tubes.
- 9 Completely fill the cell until you see solution emerge from the outlet access tube into the empty pipette tip.
- 10 Hold the outlet pipette tip vertically and make a few short strokes with the pipette plunger to gently force solution back and forth through the capillary cell. Be careful not to introduce air bubbles into the cell by letting either pipette tip go dry. This technique helps to get rid of any bubbles that might still be trapped in the capillary channel.
- 11 Make sure the reference cell is free of bubbles then remove the pipette tips and pipette from the inlet and outlet access tubes. Minimize solution spillage by covering the open end of the free pipette tip with your thumb while removing it from the access tube.
- 12 Repeat the entire process to fill the sample cell (right-most pair of access tubes).
- 13 Cap one of the access tubes with the small black access tube cap (see the figure below). Leave one access tube open on each cell to allow pressurization.



Figure 7 Cap placement on rear access tubes.

- 14 Remove excess liquid from the surrounding area outside of the access tubes using suction or a lab tissue.

Cell Pressurization and Checking for Bubbles

The last steps needed when loading the cells are to attach the pressure handle to the access flange, pressurize the cell and check for bubbles. Follow these steps, in order:

- 1 Screw the threaded ring on the pressure handle onto the access flange tightly, as shown in [Figure 8](#).



Figure 8 Attaching the pressure handle.

CAUTION: Do not pressurize the cells at this time.

MISE EN GARDE: Ne mettez pas les cellules sous pression à cet instant.

- 2 Make sure the Nano DSC is powered on.
- 3 Open the DSCRun™ program by double clicking the DSCRun icon located on the Windows desktop (if it is not already running).
- 4 Wait for the baseline to stabilize within a few tenths of a microwatt before beginning an experiment. Equilibration is complete when the μW reading (as displayed in the Output box shown in the figure below) settles to constant value.




Figure 9

- 5 Set the manostat pressure to a value of at least 3 atm (up to 6 atm for high-temperature instruments) using **Tools > Runtime variables**.
- 6 Increase the pressure by clicking on the up arrow in the DSCRun program (see the figure below) while monitoring the manostat pressure as displayed on the computer monitor.
 - If a bubble is present, the μW reading will exhibit a significant momentary change (in excess of $30\ \mu\text{W}$) as the pressure is increased from 0 to 3 atm.
 - If bubbles are indicated by an unstable μW reading as the cell pressure is applied, the cells must be refilled using the procedures outlined previously.
 - If you see little or no effect, proceed to the next step.
- 7 Set the experiment parameters and start the scans. See the online Help for more information.

An important final step in any experiment is cleaning the calorimeter. The Nano DSC is an extremely sensitive instrument that will perform well only if kept scrupulously clean. See the next section for information.

Cleaning the Nano DSC

Immediately after every experiment sequence (that is, after the baseline and sample scans are both completed for a specific experiment) the Nano DSC cells should be washed.

CAUTION: The pressure handle should not be attached to the calorimeter during cleaning.

MISE EN GARDE: Le manche à pression ne doit pas être fixé au calorimètre pendant le nettoyage.

NOTE: Do not clean the cells between the buffer/buffer and the buffer/sample scans).

The general procedure for cleaning the DSC cells with formic acid is described below.

The following cleaning process is suitable for most samples (more vigorous cleaning is accomplished by following the procedures described later in this chapter):

- 1 Wait until the calorimeter is idling at approximately 20 to 30°C, then reduce the cell pressure to zero.
- 2 Click down arrow and **remove the pressure handle** from the access flange when depressurization is complete.
- 3 Empty the cells as follows:
 - If you have a cylindrical cell, use the 1- mL syringe to remove the solutions from the cells.
 - If you have a capillary design, remove the cell outlet tube caps first and then empty the cells with the pipette or vacuum pump.

CAUTION: Ensure that the small black caps (capillary cell) and the pressure handle have been removed.

MISE EN GARDE: Assurez-vous de retirer les petits bouchons noirs (cellule capillaire) et le manche à pression.

- 4 Fill the cells with the recommended cleaning solution (50% formic acid), then run one scan cycle from 25 to 80°C at 1°C per minute. Do not install the pressure handle.
- 5 **When the temperature has returned to approximately 25°C, aspirate and discard the cleaning solution. Thoroughly rinse all areas that come into contact with corrosive chemicals. Do not let any solutions pool inside the access tube area.**
- 6 Rinse the cells with a generous volume of high quality deionized water using the appropriate equipment for your cell design. See the next two sections and follow the flushing directions for your instrument.

CAUTION: Thoroughly rinse all areas that come into contact with corrosive chemicals. DO NOT let any solutions pool inside the access tube area.

MISE EN GARDE: Rincez à fond toutes les parties qui entrent en contact avec les substances chimiques corrosives. Ne laissez jamais de mélanges de solutions à l'intérieur de la zone du tube d'accès.

Cylindrical Cell Design Flushing Equipment

After washing the cylindrical cell with formic acid, follow the instructions in this section to flush the cell.

- 1 Place the cylindrical cell cleaning adapter (shown in [Figure 10](#)) into the calorimeter's cells (as shown in [Figure 11](#)).




Figure 10 Cylindrical cell cleaning adapter and cross-section of cleaning flow.



Figure 11 Positioning the cylindrical cell flushing apparatus.

- 2 Place the inlet tubing in a large volume of water and attach the outlet tube to a vacuum pump trap flask. This allows you to flush large volumes of liquid through the cell.
- 3 If the formic acid rinse appears inadequate (as judged by noisy water/water baseline scans), follow the intense cleaning procedure as described later in this chapter.

Capillary Cell Design Flushing Equipment

After washing the capillary cell with formic acid, follow the instructions in this section to flush the cell.

- 1 Attach the Manosil silicone tubing to the access tubes for each cell in the calorimeter.
- 2 Connect one side to a large flask of water, and the other side to a vacuum pump trap flask as shown in the figure below. This allows a large volume of solution to be drawn through the calorimeter's cells.
- 3 If the formic acid rinse appears inadequate, (as judged by noisy water/water baseline scans), follow the intense cleaning procedure as described later in this chapter.



Figure 12 Capillary cell cleaning apparatus for flushing large volumes.

Software

Several software programs are available for use with the Nano DSC. The table below outlines the usage for each one. For more details on these programs, see the online help.

Table 2: Nano DSC Analysis Software

Program	Functions
DSCRun	<p>This program controls the operation of the Nano DSC and is used for data acquisition.</p> <p>The main components are:</p> <ul style="list-style-type: none">• DSC control panel• Main menu functions• Operating modes and status
NanoAnalyze	<p>Thermodynamic analysis of experimental data is conducted in this program. It allows you to calculate the partial molar and excess heat capacities of the biopolymer. Afterwards, the data can be fit to a series of models that are appropriate to the system being studied.</p> <p>The following features are present in this program:</p> <ul style="list-style-type: none">• Viewing, exporting, and importing data files.• Conversion of microwatt data into molar heat capacity. It requires a sample data file and a baseline data file for operation.• The ability to fit molar heat capacity or excess heat capacity data to one of four models. All of the models assume that there are two-state or multi-state transitions in the data. <p>Data can be exported from NanoAnalyze to Microsoft® Excel.</p>

Maintaining the Nano DSC

Maintaining the Nano DSC consists of performing balance runs when baseline shifts indicate a need, purging the instrument, cleaning the filter, lubricating the o-ring, and thoroughly cleaning the cells. This section provides information on these procedures.

Balance and Residual

NOTE: Use of the residual scan feature will not affect your data in any way.

The balance and residual scan scans serve to flatten the baselines, allowing for easier visual interpretation especially in real time or when collecting data for presentation. Perfect flatness is not a strict requirement; when baseline scans are subtracted from sample scans, the baseline shifts cancel out. Users may find it convenient to renew the baseline and residual scans from time to time, as the baseline shifts over time. This will result in water baseline data that do not fall within ± 100 μ watt range.

The measuring unit of your Nano DSC contains two cells that are closely matched in mass and internal volume. These cells are attached to the calorimeter block via the access tubes. During the manufacture of your calorimeter, these access tubes were machined to match the thermal conductivity of the sample and reference cells. Since it is impossible to achieve a perfect match of the thermal paths for the two cells, an instrument baseline, or balance run, is performed. This balance run creates a baseline for the calorimeter, which yields a heat flow signal that is near zero when both cells contain solutions with the identical heat capacity.

The balance run will give you a satisfactory baseline. However, to further flatten the baseline across the entire temperature range, you can correct the instrument baseline by running a “*residual*” scan. A residual scan is generally used for cosmetic reasons only; the real sample heat capacity data is obtained by subtracting data from the baseline, not from individual scans. Use the same scan rate and temperature limits as the experiment.

The residual feature of this software is simply a process used to adjust the linearity of the instrument baseline to display a value near zero.

The procedure for balancing the cells and creating residual scan data are as follows:

- 1 Load the cells with degassed, deionized water. Be sure you have no bubbles.
- 2 Select **Tools > Balance Cells** from the DSCRun menu bar. The calorimeter performs a 0 to 130°C scan (both up and down).
- 3 Allow the cells to equilibrate.
- 4 Select **Tools > Residual data base** from the DSCRun menu bar (see the three figures below).

Set the attributes of the **Residual data base** dialog box as shown below. (Examples only. **Rate** and **T-lower** must match the intended future experiment; **T-upper** must equal or exceed the experiment final temperature.)

- **Rate:** 1
- **T-lower:** 0
- **T-upper:** 130

- 5 Remove the previous residual values by highlighting each value on the residual data panel (**Label** box) and then clicking **Remove** (see the figures below).




Figure 13 Residual Baseline screen; Add Baseline dialog box.

NOTE: The items in the **Label** box must be removed before a new residual can be run.

NOTE: It is perfectly acceptable to run the Nano DSC with only a balance and no residual (that is, the residual values are zero (0) in the **Label** box), or a residual scan without re-balancing the instrument first.

- 6 Select **Add** (or **Rescan**) to begin. Again, the instrument will scan up and down from 0 to 130°C.
- 7 Wait for the scans to finish. Then select **File > Save** from the main menu to store the residual values.
- Your instrument is now balanced and ready for use.

A different scan rate other than 1°C/min requires a different residual scan to achieve a flat baseline. For example, if your application requires a 0.5°C/min scan rate, perform a residual scan (you do not need to run another balance scan) as described above, but change the **Rate Attribute** value from **1** to **0.5**.

Purging the Nano DSC

The Nano DSC features a sealed, nitrogen-filled canister surrounding the measurement cells in order to ensure stable operation. If baselines become unstable, it may be necessary to contact TA Instruments for service. Check for stability by monitoring the signal repeatability, especially at temperatures above 60°C. If the signal varies significantly with no pattern from one scan to the next, contact TA Instruments Service to determine if there is a need to service the instrument.

Important Note: The canister has been vacuum-purged and back-filled with dry nitrogen gas at the factory before shipment, and will only very rarely be required again in the field. Contact TA Instruments for service if the instrument does not maintain a stable signal as described above.

Cleaning O-Rings

In order to ensure a good seal for the pressure handle, you may need to occasionally remove the o-ring from the access flange with a pair of dull tweezers. Be very careful not to tear or puncture the o-ring or scratch the surface of the seating channel. Clean the o-ring groove and o-ring, then lubricate the o-ring with vacuum grease and reinstall it.

Cleaning the Instrument

Regular Cleaning

During normal operation it is recommended that you clean the Nano DSC sample and reference cells with a 50% formic acid solution between samples. Remove the sample then load the formic acid into both cells. Scan from 25 to 75 °C and back to 25 °C. After the cells have cooled to 25 °C, rinse each cell with 1 L of deionized water. See page 29 for details.

If, however, there is excessive noise in subsequent scans, or you suspect a residue or precipitate in the cells, follow the cleaning procedure below.

Thorough Cleaning

These procedures are used to more thoroughly clean the instrument if needed. Select the process that is appropriate to the suspected contaminant. Sodium hydroxide cleaning is always followed by a formic acid neutralization step.

All rinsing or flushing operations should be done with the cleaning devices described beginning on page 25. Caustic solutions should be loaded and removed with the appropriate syringe or micropipette.

Protein Deposits: Pepsin Solution

- 1 Prepare a solution containing 0.5 M NaCl, 0.1 M Acetic Acid and 1 mg/mL Pepsin.
- 2 Place this solution in the calorimeter for a minimum of three hours, preferably at 30°C. (Use fresh pepsin each time as it eats itself up in solution.)
- 3 Flush the cells with 1 to 2 L of deionized water after the solution has been in the calorimeter for three hours or more.

Mineral Deposits, Step A: NaOH Solution

- 1 Prepare a 4.0 M NaOH solution.
- 2 Place this solution in the calorimeter and run a scan from 25 to 90°C at 2°C/min. Interrupt the scan at 90°C and let it stand at this temperature overnight.
- 3 Flush the cells with 1 to 2 L of deionized water after the NaOH solution has been in the calorimeter over night.

Step B: Formic Acid Solution

- 1 Prepare a 50% formic acid solution.
- 2 Place the solution in the calorimeter and run a scan from 25 to 65°C at 2°C/min. Interrupt the scan at 75°C and let it stand at this temperature for 20 minutes.
- 3 Flush the cells with 1 to 2 L of deionized water after the 20 minutes have passed.

Grease or Oil

If you suspect grease or oil in the cells, follow these steps.

- 1 Place HPLC grade tetrahydrofuran in the calorimeter and run a scan from 25 to 50°C at 2°C/min.
- 2 Interrupt the scan at 50°C and let it stand at this temperature for 20 min.

Soapy Water Flush

- 1 Exit from the DSCRun software. (This must be done whenever you pass large volumes of solutions through the Nano DSC.)
- 2 Flush or aspirate 100 mL of soapy water (for example, SDS) through the cells.

Final Flush

Flush or aspirate 1 to 2L of deionized water through each cell.

Troubleshooting the Instrument

The following table provides tips on how to solve various problems that you may encounter when using the Nano DSC. For anything not solved below, please call TA Instruments for service.

Table 3: Common Troubleshooting

Symptoms	Possible Problem	Action/Remedy
Random sharp spikes in scan traces.	<ul style="list-style-type: none"> Air bubbles in sample. Air bubbles in Nano DSC. 	<ul style="list-style-type: none"> Degas sample. Fill cells correctly so no air pockets develop.
Noisy scans.	<ul style="list-style-type: none"> Contamination from previous sample. Precipitation of sample. Sample contamination. Canister may need to be evacuated. 	<ul style="list-style-type: none"> Clean cells thoroughly. Clean cells with 4 M NaOH solution. Rinse thoroughly. Avoid running a sample that precipitates or is a suspension. Check sample for impurities. Check syringe and needle for contaminants. This service is rarely needed. Consult with TA Instruments representative before proceeding.
Scans have a sharp jagged rise between 90 and 100°C.	Sample is boiling out of cells.	Cells may not be pressurized. Check for possible causes: <ul style="list-style-type: none"> Pressure handle not installed Manostat not adjusted to elevated pressure Missing, torn, or contaminated pressure handle o-ring
Only the last scan is displayed during the run.	View option setting.	Select All scans under the View option on the toolbar.
Water-water scans do not fall within the $\pm 100 \mu\text{W}$ range.	<ul style="list-style-type: none"> Balance has degraded over time. The cells are unevenly filled. 	<ul style="list-style-type: none"> Rebalance the Nano DSC. Run a residual. Refill the Nano DSC carefully with the same amount of solution in each cell.
Consecutive base-lines shift by more than $3 \mu\text{W}$.	<ul style="list-style-type: none"> Cells require cleaning cure. Canister may need to be evacuated. Pressure leak. 	<ul style="list-style-type: none"> Select and perform appropriate cleaning procedure. This service is rarely needed. Consult with TA Instruments representative before proceeding. Check the pressure handle and o-ring.
Nano DSC does not hold constant pressure.	<ul style="list-style-type: none"> O-ring is dry, contaminated, or damaged. Pressure handle is not screwed on tightly. 	<ul style="list-style-type: none"> Carefully remove the o-ring (do not use sharp metal tools). Inspect, wipe clean, and lubricate the access o-ring. Replace if necessary. Tighten the pressure handle.

Table 3: Common Troubleshooting (Continued)

Symptoms	Possible Problem	Action/Remedy
Nano DSC does not pressurize	<ul style="list-style-type: none">• No power to instrument (Green light on front panel should be illuminated)• Pressure motor piston not in initial position	<ul style="list-style-type: none">• Check power cord, surge suppressor strip, and electrical service.• Remove pressure handle and reset by clicking on depressurize control (“down” arrow)

Appendix: A

Pressure Perturbation Experiments

Colette Quinn, Ph.D. and Lee Hansen, Ph.D.
TA Instruments, 890 West 410 North, Lindon, UT 84042, USA

Introduction

Pressure Perturbation Calorimetry (PPC) is a technique for determination of thermal and volumetric properties of macromolecules. These properties include the partial molar thermal cubic expansion coefficient ($\bar{\alpha}$, Equation 1, [1]) and the relative partial molar volume ($\Delta\bar{V}/\bar{V}$, Equation 2, [2]) of the macromolecule as functions of temperature. These parameters are associated with hydration of the macromolecule [3, 4]. The method consists of measuring the heat effect, ΔQ , from applying a pressure change, Δp , to a solution of the macromolecule. $\bar{\alpha}$ and $\Delta\bar{V}/\bar{V}$ are then calculated by the following equations.

$$\bar{\alpha} = \left(\frac{1}{\Delta V}\right) \left(\frac{\partial V}{\partial T}\right)_p = \frac{\Delta Q}{VT\Delta p} \quad \text{Equation 1}$$

$$\frac{\Delta\bar{V}}{\bar{V}} = \int_{T_o}^{T_e} \bar{\alpha} dT \quad \text{Equation 2}$$

The Nano DSC is currently equipped with the necessary components to perform PPC without further modification. PPC is typically performed on a DSC that is capable of increasing and decreasing the pressure on a sample during a DSC temperature scan or during isothermal measurements. During a temperature scan, the pressure change must occur faster than the heating or cooling of the system (3, 4). This note is a short tutorial on the required steps to set-up an experiment on the Nano DSC with DSCRun software and then on analysis of the resulting PPC data with NanoAnalyze software to evaluate $\bar{\alpha}$, which denoted as TEC in the software.

PPC with DSCRun

Turn on the Nano DSC, then open DSCRun and allow the system to stabilize. In the DSCRun software open **set runtime variable** under the **Tools** menu.

On the **Experimental Method** tab, check either scanning or isothermal and enter the temperature limits and scan rate; see [Figure 14](#). If scanning, a slow rate is required because the calculations assume isothermal conditions; see the definition of $\bar{\alpha}$ in Equation 1 (4). If the scan rate is slower than the instrument feedback then it can be safely assumed isothermal during the pressure steps. A rate of 0.1 °C/min is recommended. This scan rate will typically satisfy the isothermal requirement without scanning so slowly that resolution of the heat rate is unduly sacrificed.

Choose the pressure function to be applied during the experiment. See [Figure 14](#).




Figure 14 Screen shot of DSCRun highlighting the different pressure parameters available.

Three different pressure control options are available: linear, sine, or ramp (labeled **jump** in the software). The jump pressure function is usually used. Once the desired pressure function is selected and all **Temperature Parameters** are set for the scan, click **Add To Experiment Method**. This opens a dialog box that enables setting the **Jump Pressure Function** (see [Figure 15](#)).




Figure 15 Screen shot of dialog box that opens after a pressure function has been selected. (A indicates atmospheres and s indicates seconds.)

Note that in scanning experiments, pressure can be designated as a function of time or temperature. The term "step" in this box refers to the pressure step. Although the instrument is capable of operating from 0 to 6 atm, relative to ambient pressure, the recommended pressure change is from 1 to 5 atm with a ramp/step time of 600 s (4). The minimum pressure must exceed a pressure that allows boiling at any temperature during the scan. For most experiments, the default motor speed of 28×10^{-1} ms/step is sufficiently slow to satisfy the requirement of isobaric conditions in Equation 1.

Experimental Set-up

A water-water scan can serve as a test experiment to assess any asymmetry between the reference and sample cells. For this test, water is loaded into both the reference and sample cell. The pressure is stepped from 1 to 5 atm, starting at low pressure with time steps of 600 s while scanning at 0.1 °C/min from 20 to 80 °C. During the experiment pressure is applied equally to both sample and reference cells, so any deviation of ΔQ from zero is solely due to any asymmetry between the cells.

For typical solutions of macromolecules, a concentration of 10 mg/mL should generate good data. Due to the slow scan rate, concentration requirements for PPC experiments are greater than for a normal DSC scan. The PPC experiment consists of two parts, a background run and a sample run. Background data with buffer in both reference and sample cells should be collected under identical PPC conditions as the sample. Sample data must be collected with the same buffer in the sample cell. The basis for this method and the calculations done in NanoAnalyze are as follows (Equations 3-5).

$$\frac{\Delta Q}{VT\Delta p} = \alpha_{measured} = (\alpha_{sample} - \alpha_{reference}) + \alpha_{asymmetry} \quad \text{Equation 3}$$

$$\alpha_{asymmetry} = (\alpha_{buffer \text{ in sample cell}} - \alpha_{buffer \text{ in reference cell}}) \quad \text{Equation 4}$$

$$\bar{\alpha}_{macromolecule} = (\alpha_{sample} - \alpha_{reference}) = \alpha_{measured} - \alpha_{asymmetry} \quad \text{Equation 5}$$

$\alpha_{measured}$ is calculated from the data with sample solution in the sample cell and buffer in the reference cell and $\alpha_{asymmetry}$ is calculated from the data with buffer in both cells.

PPC Fitting with NanoAnalyze

The data files from the background PPC run and the sample PPC run can be opened in NanoAnalyze without modification. When the raw data files are initially opened, a dialogue box will open with a prompt: **Analyze pressure data for Experiment Step 1?** (Figure 16). In this context, "step" indicates a temperature scan or isothermal run. This request will be made for each temperature scan or isothermal run if the box at the bottom is not checked.




Figure 16 Dialogue box prompt when opening PPC data with NanoAnalyze.

When the files are opened, extra data analysis tabs are automatically made accessible ([Figure 17](#)).




Figure 17 Screen shot highlighting the extra tabs accessible after opening PPC data on NanoAnalyze.

In the **Input Data** tab, several key values must be entered for TEC calculations, i.e., sample cell volume, concentration, and molecular weight in kilodaltons. The data file is loaded with the **Add file** icon. The background data is dragged and dropped into the **drop baseline here** box ([Figure 17](#)). MHC (molar heat capacity), the second tab, is not used in PPC calculations. The third tab from the left, **Integration Baseline**, gives the integrated area under each pressure events. The baseline for integration of the heat effects (ΔQ) from the pressure changes can be either a straight-line (default) or a sigmoid (check **Sigmoid Baseline** box) ([Figure 18](#)).




Figure 18 PPC sample and background data in the Integration Baseline window.

In the **Area graphs** tab (Figure 19), background and sample areas are plotted. The baseline area data can be fit to a polynomial, labeled blank fit order in the software (as is the case with TEC Baseline and TEC tabs).




Figure 19 Area graphs window for a sample.

The **TEC Baseline** is used for calculating α for pure liquids, and is not used for calculation of $\bar{\alpha}$ for solutes.

The **TEC** tab plots the thermal expansion coefficients in the top window (Figure 7) and the integrated area in the bottom window. This integrated area is the relative change in volume ($\Delta\bar{V}/\bar{V}$), according to the relationship of Equation 2 (3).




Figure 20 The TEC (α) for a sample.

Numerical values for all of the thermal expansion data are displayed in the **Table** tab and can be copied to a clipboard and imported to another program if further manipulation is desired.

References

- 1 Randzio, S.L. "Comments on "volumetric studies of aqueous polymer solutions using pressure perturbation calorimetry . . . " [Macromolecules 34 (2001) 4130]" *Thermochimica Acta* 2003 398, 75-80.
- 2 Mitra, L.; Smolin, N.; Ravindra, R.; Royer, C.; Winter, R. "Pressure perturbation calorimetric studies of the solvation properties and the thermal unfolding of proteins in solution-experiments and theoretical interpretation" *Phys. Chem. Chem. Phys.* 2006 8, 1249-1265.
- 3 Rosgen, J. "Pressure-Modulated Differential Scanning Calorimetry: Theoretical Background" *Anal. Chem.* 2006 78, 991-996.
- 4 Dragan, A.I.; Russell, D.J. Privalov, P.L. "DNA Hydration Studied by Pressure Perturbation Scanning Microcalorimetry" *Biopolymers* 2008 91, 95-101.

Appendix: B

Guardian Option

Guardian Option

The Notifications tab is where Guardian option features reside. Guardian is an implementation of 21 CFR Part 11 regulations. Detailed information about Guardian is available in the document “Guardian for Microcalorimetry Software.”




Figure 21 Guardian tab.